


EXPERIENCE THE **HYBRID**

FLEXIBLE. DURABLE. POWERFUL

The new PowerFlex **Hybrid*** High Power Holmium fiber is a breakthrough in fiber technology. Its innovative design provides the superior flexibility of a 200 μ fiber while maintaining the stone treatment benefits of a larger 272 μ round tip fiber. That combined with the new and improved internal coating facilitates safe energy transmission under a 1 cm bend using higher power settings.¹



13%
greater
surface area
than a regular
200 μ fiber

POWERFlex™-HYBRID

HOLMIUM ROUND TIP LASER FIBER


CLARION[®]
MEDICAL TECHNOLOGIES

800.668.5236
www.clarionmedical.com

ROUND TIP

Designed to glide freely through tightly deflected scopes.

ARMOURTEC COATING

An advanced coating for a durable design allowing the fiber to withstand a 1.0 cm bend under high power settings.



POWERFlex™-HYBRID

HOLMIUM ROUND TIP LASER FIBER

ZIRCONIA CERAMIC CONNECTOR

Exclusive ceramic connector safely diffuses errant energy to provide a high level of reliability and eliminates damage to the debris shield and laser optics.

HYBRID TIP

Larger more durable surface area for faster and more effective clinical results.

UNIQUE PACKAGING

Patented dispensing tray protects fibers against pinch points and allows for safe and controlled removal.

POWERFLEX HYBRID ORDERING INFORMATION

| Fiber Size | Fiber Size at Distal Tip | Description | Maximum Outer Diameter | Units per Box | Part # |
|------------|--------------------------|--|------------------------|---------------|----------------|
| 242µ | 272µ | PowerFlex Hybrid 200 Single-Use High Power Holmium Laser Fiber | 450µ | 5 | PF200SE-RT-5PK |

* Patent Pending
1. Bench Testing on File